

## QuantIFERON-TB Gold Plus: Results and Interpretation

Nil	TB1 - Nil	TB2 - Nil	Mitogen - Nil	Qualitative Result	Interpretation
≤8.0 IU/mL	≥0.35 and ≥25% of Nil	Any	Any	Positive	<i>M. tuberculosis</i> infection likely
	Any	≥0.35 IU/mL and ≥25% of Nil			
	<0.35 IU/mL OR ≥0.35 IU/mL and <25% of Nil		≥0.5 IU/mL	Negative	<i>M. tuberculosis</i> infection NOT likely
>8.0 IU/mL	Any			Indeterminate	Likelihood of <i>M. tuberculosis</i> infection cannot be determined

QFT-Plus is considered positive if either one or both Tb Ag tubes are positive

Interpretation table adapted from Mayo Medical Laboratories

(<https://news.mayomedicallaboratories.com/2018/01/22/quantiferon-tb-gold-plus-hot-topic/>)