

Breo Ellipta



Breo Ellipta® (fluticasone furoate & vilanterol)

100/25mcg (blue label)

200/25mcg (yellow label)



Medicine	Need a spacer? spacer?	My dose
Controller (dry powder inhaler)	No	Check Asthma Action Plan

1. Holding Ellipta upright, slide the cover down to expose the mouthpiece and to hear a click
2. Exhale away from device
3. Place mouth over mouthpiece with a tight seal (do not block vent with fingers)
4. Take one long steady deep breath in through your mouth.
5. Hold breath as long as you can, about 3-4 seconds
6. Remove Ellipta from mouth and exhale slowly away from device.
7. Close cover back over mouthpiece as far as it will go
8. Repeat steps 1-7 for each dose



Use the camera on your phone to scan this QR code for the Breo Ellipta Website:

1/6/2021

There is information about:

- ◆ How to use the Ellipta
 - ◆ Click on “Patient Information”
 - ◆ Click on “See a demo on how to use Breo Ellipta”